



Important hints

This investment must not come into contact with plaster or investments containing plaster. Use separate mixing tools.

If the wax model is treated with a surface tension reducing wetting agent, let it dry completely before investing. Never use compressed air for drying.

Attention

This investment contains quartz and cristobalite, therefore avoid inhaling dust when handling the powder

Packing units

Portable cartons with foil insert	20 kg
Pouches	180/400 g

Storage

In well closed, vapour-proof packages. Pay attention to the expiry-date on the packages.

The above recommendations are given to the best of our knowledge. We grant the quality of our products according to our specification. Any further liability cannot be accepted since the proper application of our products is outside of our control.